Lab 5 Reverse Transcription Polymerase Chain Reaction (RT-PCR)



30 - 40 cycles of 3 steps :



5' TTTTTTTTT

5'

3"

3'

Step 1 : denaturation

1 minut 94 °C

Step 2 : annealing

45 seconds 54 °C

forward and reverse primers !!!



ST

3'

You will need : -RNA sample----- mRNA-----cDNA - primers -dNTP - buffers RNase and RNase inhibitor -RT enzymes

RT PCR: Reverse Transcription Polymerase Chain Reaction



One step and two step prosedures

-One step :RT and PCR performed consecutively In a single tube

-Two step: RT and PCR performed in seperate tubes

Advantages of one step procedure: -Minimize time requires -Reduces risk of contamination -Improve sensitive and specificity

Advantages of two step procedure: -Allows optimal reaction conditions -Provide maximum flexibility -Amplifies long sequences **Choosing the RT-PCR enzymes:** -Maximum lenght of tamplet that can be transcribed into full-lenght cDNA -Temperature optimum -Snsetivity -Specificity -Rnase H activity

Choosing the primers for TR: -Oligo (dt)n -Anchored oligo (dt)n -Random hexamers -Sequence-specific

RT- PCR

Using Reverse Transcriptase to convert mRNA into complementary DNA (cDNA) which was then amplified by PCR and, again analyzed by agarose gel electrophoresis.

Reverse Transcriptase-PCR analysis of mRNA is often referred to "RT-PCR" which confused with "real time-PCR".

Complementary DNA

cDNA is synthesized from mature mRNA using the enzyme reverse transcriptase and the enzyme DNA polymerase. This enzyme operates on a single strand of mRNA, generating its complementary DNA based on the pairing of RNA base pairs (A, U, G and C) to their DNA complements (T, A, C and G respectively).



Reverse transcriptase uses a single-stranded RNA template to create a double-stranded cDNA.

DNA library

is a collection of cloned DNA fragments. <u>There are</u> <u>two types of DNA library:</u>

- <u>Genomic library</u> contains DNA fragments representing the entire genome of an organism.

 <u>cDNA library</u> contains only complementary DNA molecules synthesized from mRNA molecules in a cell.
It contain all the coding information (the genetic code that produces proteins) from a given cell population.

cDNA Library Technique:

1- After the cDNA is synthesized, it is cloned into expression vectors or plasmids.

2- These plasmids each containing one cDNA are transformed into bacterial competent cells .

3- These plasmids are amplified in the growing bacteria.

4- The bacteria clones are then selected so that only bacteria containing the plasmid will survive.

5- This is commonly done through antibiotic resistance selection.

6- Once the bacteria are selected, stocks of the bacteria are created which can later be grown and sequenced to compile a cDNA library.

