
Lab.343

Experiment (9)

Molecular weight of a polymer from viscosity measurements


polymers

are a large class of materials consisting of many small molecules (called *monomers*) that can be linked together to form long chains.




Types of polymers:


- 1-Natural
 - ex. Natural rubber, <u>DNA</u> and <u>proteins</u>
- 2-Industrial(synthetic)
 - ex. <u>plastics</u> such as <u>polystyrene</u>



Polystyrene > $(C_8H_8)_n$



Polymerization











The ratio of the viscosity (n) of a solution polymer molecules to the viscosity (n_o) of the solvent related to the molecular weight M of the polymer

$$\frac{(\eta / \eta_{o} - 1)}{C} = K M^{\alpha}$$


K: is a constant for any given type of polymer.
a: is a function of the geometry of the molecule.
C: is the number of grams of polymer in 100 ml of solution.

(η / $\eta_{\circ}\,$ - 1) is known as the specific viscosity

 $\eta_{sp}/C = K M^{\alpha}$

- This equation is only valid for very dilute solution.
- The graph which is drawn of η_{sp}/C vs C.
- Intrinsic viscosity [n]

 $[\eta] = \lim_{c \to 0} \eta_{sp}/C$



$$[\eta] = K M^{\alpha}$$

Procedure:


- Prepare the stock solution of the polystyrene in toluene by weighting 2.5gm of the polystyrene dissolved in 200ml toluene.
- 2. Make 5 different concentration by diluting the stock solution to 2,4,8 and 16 times.
- 3. Measure the average time of flow (t) for each of the above 5 concentrations using viscomter at $25^{\circ}C$.
- 4. Measure the average flow time of toluene (t $_{\circ}$) at 25°C.

Results

С	t=	t。	$\eta/\eta_{o} = t/t_{o}$	$\eta_{sp} = \eta/\eta_{o} - 1$	η_{sp}/C
	η	= η _°			
0.5					
.25					
0125					
.062					
.03125					

$[\eta] = K M^{\alpha}$

Polymer-solvent system	K	α
Poly styrene-Toluene	1.05 x 10 ⁻⁴	0.725





.....=K M ^{0.725} M ^{0.725} =/ K

 $0.725 \ln M = \ln (..../K)$

ln M = /0.725M= g/mol

**Plot { 1/C ln (η / ηo) } vs. C. and compare the results