

Lab.343



Experiment (3)

**Determination of the saponification rate
constant of ethyl acetate in alkaline
medium**

A. PURPOSES

The purpose of this experiment is to show the reaction of ethyl acetate saponification by hydroxide ions :



is a second-order reaction. In addition, it is also determined the reaction rate constants. This determination is done by titration method.

BASIC THEORY

To determine the rate of a given chemical reaction, it should be determined how fast changes occur in the concentration of reactants or products.

In general, if a reaction occurs $A \rightarrow B$, the initial substance A and substance B did not exist.

After some time, the concentration of B will increase while the concentration of A will decrease.

Rate law can be determined by conducting a series of systematic experiments on the reaction $A + B \rightarrow C$, to determine the reaction order with respect to A the concentration of A is fixed while B concentration varied then determined the rate of the reaction on the concentration variation. As for determining the order of the reaction B, the concentration of B is fixed while the concentration of A varied then measured the rate of the reaction on the concentration variation



Hydrolysis is a chemical decomposition involving breaking of a bond and the addition of elements of water

In this hydrolysis of ester (ethyl acetate) with an alkali (sodium hydroxide), HCl was used as catalyst to accelerate it.



t= 0	a	b	0	0
t= t	a-x	b-x	x	x

$$\begin{aligned} R = dx/dt &= K [\text{CH}_3\text{COOC}_2\text{H}_5][\text{NaOH}] \\ &= K [a-x][b-x] \end{aligned}$$

Second order reaction equation

If

• $[a] \neq [b]$

Or

• $[a] = [b]$

• $[a] \neq [b]$

$$dx/dt = K (a-x) (b-x)$$

$$dx / (a-x) (b-x) = K dt \quad \textit{by integration}$$

$$K = [2.303 / t (a-b)] [\log b(a-x) / a(b-x)]$$

• $[a] = [b]$

$$dx/dt = K (a-x) (b-x)$$

$$dx/dt = K (a-x)^2$$

$$dx / (a-x)^2 = K dt$$

by integration : $\int_0^x dx / (a-x)^2 = \int_0^t K dt$

$$1/a-x = Kt + \text{constant}$$

At $t=0$, $x=0$  constant = $1/a$

$$1/a-x = Kt + 1/a$$

At $t = t_{1/2}$, $x = a/2$

$$t_{1/2} = 1/ak$$

unit of $t_{1/2}$ = time

$$K = 1/a t_{1/2}$$


unit of K = $L \text{ mol}^{-1} \text{ time}^{-1}$

Prepare:

concentration	Volume	Liquid	Solid
0.05M	100ml	$\text{CH}_3\text{COOC}_2\text{H}_5$	NaOH
0.025M	250ml	HCl	NaOH

Procedure:


1-




**Take room temperature
10 minutes**

2-pour base to ester (fast) then open stop watch with strong shaking.


3-At time = 2,4,7,... take 10 ml of mixture



4-



time	$V(\text{NaOH}) \equiv X \equiv \text{HCl}$	$a-x = 10 - x$	$1/a-x$
2			
4			
7			
10			
14			
18			
23			
28			
38			



➤ $\text{HCl react} \equiv \text{free NaOH (unreacted)} (a-x)$

➤ $\text{HCl unreacted} \equiv V \text{ NaOH} \equiv \text{reacted NaOH} + \text{reacted ethyle acetate} \longrightarrow \text{product} \equiv X$