

General Biochemistry

Organized by

Lecturer\ **Duaa Zahim**

Lecturer\ **Dalal Al-Saedi**



Name:

Computer no.

Section

Wear

Personal
Protective
Equipment
(PPE)



Avoid

Eating and Drinking



Messing



Change dropper from
bottle to another

Keep

Lab notebook



Calm



Lab safety and
follow instruction



State of Materials

Solid

Liquid

In the laboratory

Terms

Concentrated

Diluted

lab Equipments

Test Tubes



Test Tube Holder



Test Tube Rack



Test Tube Brush



Dropper



Bottels



Microscope Slides



Chemicals

Acids

Bases

Reagents

lab Instruments

Water Bath



Light Microscope



Biosafety cabinet



Vortex



Hazards



Flammable



Harmful



Corrosive



Toxic



Explosion



Oxidizer

Table of content

| # | | Experiment |
|----|----------------------|---|
| 1 | Carbohydrates | Qualitative analysis of carbohydrates (monosaccharides) |
| 2 | | Qualitative analysis of carbohydrates (disaccharides and polysaccharides) |
| 3 | | Scheme of carbohydrates |
| 4 | | First Exam |
| 5 | Lipids | Qualitative Determination of Lipids |
| 6 | | Scheme of lipids |
| 7 | Protein | Qualitative Determination of proteins |
| 8 | | Scheme of proteins |
| 9 | | Scheme of unknowns |
| 10 | | Second Exam |
| 11 | | Final Exam |

Qualitative analysis of carbohydrates



Introduction of carbohydrates

Carbohydrates are organic molecules which contain carbon, hydrogen and oxygen and are the most abundant class of organic compounds found in living organisms.

Carbohydrates are produced from carbon dioxide and water by plants through the process of **photosynthesis**. They are easily digested by animals where they are converted back into carbon dioxide and water, with a concurrent release of energy. The formulas of many carbohydrates can be written as carbon hydrates, $C_n(H_2O)_n$, hence their name.




The function of carbohydrates

It is the **main energy source** for animals as they breakdown carbohydrates during the process of metabolism to release energy. Carbohydrates also serve as a **structural material (cellulose)**, a component of the **energy transport** compound **ATP**, recognition sites on cell surfaces, and one of three essential components of DNA and RNA.



The classification of carbohydrates

1- **Simple carbohydrates**, often called **monosaccharides**, contain one saccharide unit and cannot be broken down into smaller carbohydrates. These are the simplest form of sugar and are usually colorless, water soluble and crystalline solids, e.g. **glucose**, **fructose**, **galactose** etc.



The two descriptors are commonly combined into a single term like "aldohexose" for an aldose or "ketohehexose" for ketose that also is a hexose. See below some examples of monosaccharides.




Figure 1: The classified of monosaccharides according to the number of carbon atoms and functional groups.

2- Disaccharides contain **two monosaccharide** unit and **can be** broken into two monosaccharide units by hydrolysis and lost one molecule of water. See below some examples of disaccharides, e.g. **maltose**, **lactose** and **sucrose**.




Figure 2: The example of disaccharides

3- **Oligosaccharides** contain 3-10 monosaccharide unit and can be broken into 3-6 monosaccharide units by hydrolysis and lost one or more molecule of water. See below some examples of oligosaccharides, e.g. **Raffinose**.




Figure 3: The example of oligosaccharides

4- **Polysaccharides** can contain over 10 or more monosaccharide unit and can be broken into 7 or more monosaccharide units by hydrolysis and lost one or more molecule of water. See below some examples of polysaccharides, e.g. **starch**.




Figure 4: The example of polysaccharides



Lab experiments

A- Physical properties

Solubility- color- shape.

B- Chemical properties



Principle

Presence of carbohydrates can be confirmed **qualitatively** by several tests. Due to the presence of different numbers of sugar units specific carbohydrates exhibit typical color reactions that form the basis for their identification. In the following tests will be performed for the qualitative analysis of carbohydrates:

1- The Molisch test

General test for **all carbohydrates**, Monosaccharides give a rapid positive test. Disaccharides and polysaccharides react slower.

Reactions:

The test reagent dehydrates **pentoses** to form **furfural** (first reaction) and **dehydrates hexoses** to form **5-hydroxymethyl furfural** (second reaction). The furfurals further react with **α -naphthol** present in the test reagent to produce **a purple product**





Figure 5: The reaction of molisch test

How to perform the test

- 1- Two ml of a sample solution is **placed** in a test tube.
- 2- Two drops of the molisch reagent (a solution of α -naphthol in 95% ethanol) is **added**.
- 3- The solution is then **poured** slowly into a tube containing two ml of **concentrated** sulfuric acid so that two layers form.



Result

Positive test is indicated by the formation of a **violet ring** at the interface of two layers.



Negative test



Positive test

2- The Benedict test

It is used to **differentiate** between the reducing and non-reducing sugar.

Reactions

Benedict's test is **based** upon the participation of the **aldehyde** and **ketone groups** in a chemical reaction.

Reducing sugars are **oxidized** by the copper ion in **alkaline medium** to form a carboxylic acid and a reddish precipitate of copper (I) oxide.




Figure 5: The reaction of benedict's test

How to perform the test

- 1- **One ml** of a sample solution is **placed** in a test tube.
- 2- **Two drops** of the Benedict's reagent (a solution of sodium citrate and sodium carbonate mixed with a solution of copper sulfate) is **added**.
- 3- The solution is then **heated** in a boiling water bath for **two minutes**.



Result

Show **positive test** for **reducing sugars**, is indicated by the formation of an **orange-brown precipitate** within two minutes.



Negative test



Positive test

3- The Fehling's test

It is also used to **differentiate** between the reducing and non-reducing sugar.

Reactions:

Reducing sugars are **oxidized** by the copper ion in **alkaline medium** to form a carboxylic acid and a reddish precipitate of copper (I) oxide. Fehling's reagent is commonly used for reducing sugars but is known to be not specific for aldehydes.




Figure 6: The reaction of Fehling test

How to perform the test:

- 1- **One ml** of a sample solution is **placed** in a test tube.
- 2- **One drops** of the fehling's reagent (fehling's A, uses 7 gm $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$ dissolved in distilled water containing 2 drops of dilute sulfuric acid. Fehling's B, uses 35 gm of potassium tartarate and 12 gm of NaOH in 100 ml of distilled water) is **added**.
- 3- The solution is then **heated** in a boiling water bath for **two minutes**.



Result:

Show positive test for **reducing sugars**, is indicated by the **formation of a red- brown precipitate** within two minutes.



Negative test



Positive test



Notes

Sucrose is **non-reducing sugar** (does not react with Fehling's or Benedict's reagent) because the anomeric carbon of glucose is involved in the glucose-fructose bond and hence is not free to form the aldehyde in solution. Although the glucose and fructose alone are **reducing sugar**.

4- The Barfoed's test

The Barfoed's test is used to **distinguish** between reducing monosaccharides and disaccharides.

Reactions

Reducing monosaccharides are **oxidized** by the **copper ion** in acidic solution to form a carboxylic acid and a **reddish precipitate** of copper (I), oxide **within three minutes**. **Reducing disaccharides** undergo the same reaction, but do so at **a slower rate**.




Figure 7: The reaction of Barfoed's test

How to perform the test

- 1- **One ml** of a sample solution is **placed** in a test tube.
- 2- **Two ml** of the Barfoed's reagent (a solution of cupric acetate and acetic acid) is **added**.
- 3- The solution is then **heated** in a boiling water bath for **1-3 minutes**.



Result

Show positive test for **reducing monosaccharides**, is indicated by the formation of a **reddish precipitate** in the bottom of tube within 1-3 minutes.



Negative test



Positive test



Notes

In the Barfoed's test:

The **water bath** must be **boiling**.

The **sugar solution** must be **concentrated**.

The **time taken** to formation of reddish precipitate in the bottom of tube must be **calculate carefully**.

5- The Seliwanoff's test

The Seliwanoff's test is used to **distinguish** between **ketoses** and **aldoses**. If a sugar contains a ketone group it is called ketose and if it contains an aldehyde group then it is called aldose.

Reactions

The test reagent **dehydrates** ketohexoses to form 5-hydroxymethyl furfural. 5-Hydroxymethyl furfural further **reacts** with resorcinol present in the test reagent to produce a **red product** within five minutes. Aldohexoses react to form the same product, but do so more slowly.




Figure 8: The reaction of Seliwanoff's test.

How to perform the test

- 1- One half ml of a sample solution is **placed** in a test tube.
- 2- Two ml of the Seliwanoff's reagent (a solution of resorcinol and HCl) is **added**.
- 3- The solution is then **heated** in a boiling water bath for **five minutes**.

The **sugar solution** in Seliwanoff's must be **concentrated**



Result

Show **positive test** for **ketoses**, is indicated by the **formation of a red color**.



Negative test



Positive test

6- The Hydrolysis test for sucrose

Sucrose is the only **non-reducing sugar** so it does not reduce the alkaline Cu solutions (Fehling and Benedict). The sucrose first must be hydrolyzed to its component and then test it.

How to perform the test

- 1- **Two ml** of a sample solution is **placed** in a test tube.
- 2- **Two drops** of concentrated hydrochloric acid HCl is **added**.
- 3- The solution is then **heated** in a boiling water bath for **15 minutes**.



Result

Show positive test for **Sucrose**, is indicated by the formation of a light orange color.



Negative test



Positive test

7- The Iodine/Potassium Iodine test:

Iodine forms colored complexes with polysaccharides. The color of the complex depends upon the three dimensional structure of the polysaccharide. **Starch** is a coiled structure which turns **blue** when bound to Iodine.

How to perform the test:

- 1- **Two ml** of a sample solution is **placed** in a test tube.
- 2- **Two drops** of Iodine/Potassium Iodine solution and **one ml** of water are **added**.



Result

Show **positive test** for **Starch**, is indicated by the **formation of a blue- black color**.



Negative test



Positive test

8- The hydrolysis test for starch

The experiment illustrates the **conversion** of **starch** to a **reducing sugar (maltose)** by the action of sulfuric acid at **boiling point**. The longer the starch is exposed to the acid the **further** hydrolysis proceeds.

How to perform the test

- 1- **One ml** of starch solution + **One ml** of sulfuric acid are **placed** in a test tube.
- 2- The solution is **heated** in a boiling water bath for **15 minutes**.
- 3- **Add drop** from (Iodine/ Potassium Iodine) \longrightarrow Solution with **yellow color**.
- 4- **One ml** of yellow solution + **One ml** of Benedict or fehling reagent, then **heat** in boiling water bath for 1 minutes.



Result

Show positive test for **Starch**, is indicated by the formation of a **yellow color**, and then converted to **orange precipitate**.

9- The Osazone test

The Phenylhydrazine ($C_6H_5NHNH_2$) **reacts** with carbons no. 1 and no. 2 of reducing sugars to form derivatives called **osazones**. The formation of these distinctive crystalline derivatives is useful for comparing the structures of **monosaccharides** and **disaccharides**.

Reactions

The osazone reaction was developed and used by Emil Fischer to **identify** aldose sugars differing in configuration only at **the alpha-carbon**. The upper equation shows the general form of the osazone reaction, which effects an alpha-carbon **oxidation** with formation of a bis-phenyl hydrazone, known as an osazone. Application of the osazone reaction to **D-glucose** and **D-mannose** demonstrates that these compounds differ in configuration only at C-2.




Figure 9: The reaction of Osazone test

How to perform the test

- 1- One gm of the sugar sample is **placed** in a test tube.
- 2- Two gm of phenylhydrazine hydrochloride + 3 gm of crystallized sodium acetate, and 4 mL of distilled water are **added**.
- 3- **Place** the test tube in a beaker of boiling water bath.
- 4- Note the time that the test tube was immersed and the time of the precipitation.
- 5- After 20 -45 min, **remove** the test tube from the hot water bath and set it aside to cool.
- 6- A small amount of the liquid and solid is **poured** on a watch glass.
- 7- Tip the watch glass from side to side to spread out the crystals, and **absorb** some of the mother liquid with a piece of filter paper, **taking care not to crush or break up the clumps of crystals**.
- 8- **Examine** the crystals under a low-power microscope (about 80-100×), and **compare** with photomicrographs.



Result

Show positive test for **all carbohydrates** except sucrose and starch.

The formation of precipitate for **monosaccharides** (glucose or fructose) **before cooling**, while the formation of precipitate **after cooling** for **disaccharides** (maltose or lactose).



Glucosazone + Fructosazone
(Needle shaped)



Maltosazone
(Sunflower shaped)



Lactosazone
(Cotton shaped)



RESULTS & LAB REPORT

Unknown (glucose)

Physical properties

Chemical properties

| Test | Observation | Result |
|------|-------------|--------|
| | | |



RESULTS & LAB REPORT

Unknown (fructose)

Physical properties

Chemical properties

| Test | Observation | Result |
|------|-------------|--------|
| | | |



RESULTS & LAB REPORT

Unknown (sucrose)

Physical properties

Chemical properties

| Test | Observation | Result |
|------|-------------|--------|
| | | |



RESULTS & LAB REPORT

Unknown (lactose)

Physical properties

Chemical properties

| Test | Observation | Result |
|------|-------------|--------|
| | | |



RESULTS & LAB REPORT

Unknown (maltose)

Physical properties

Chemical properties

| Test | Observation | Result |
|------|-------------|--------|
| | | |



RESULTS & LAB REPORT

Unknown (starch)

Physical properties

Chemical properties

| Test | Observation | Result |
|------|-------------|--------|
| | | |



RESULTS & LAB REPORT FOR SCHEME

Unknown (1)

Physical properties

Chemical properties

| Test | Observation | Result |
|------|-------------|--------|
| | | |



RESULTS & LAB REPORT FOR SCHEME

Unknown (2)

Physical properties

Chemical properties

| Test | Observation | Result |
|------|-------------|--------|
| | | |



RESULTS & LAB REPORT FOR SCHEME

Unknown (3)

Physical properties

Chemical properties

| Test | Observation | Result |
|------|-------------|--------|
| | | |

Qualitative Determination of Lipids

Introduction

The term “**lipids**” applies to a class of compounds that are **soluble** in nonpolar organic solvents (e.g. alcohol, benzene, ether and chloroform) and **insoluble** in polar water. Lipids contain carbon, hydrogen and oxygen but have far less oxygen proportionally than carbohydrates.

Lipids are an important part of living cells. Together with carbohydrates and proteins, lipids are the main constituents of plant and animal cells.

Fatty acids: They are the building blocks of lipid. Fatty acid have a **long hydrocarbon chain** containing a carboxyl group at the end. They are divided into: saturated fatty acids and unsaturated fatty acids (unsaturated contain double bonds). The general formula for fatty acids { $\text{CH}_3 (\text{CH}_2)_n \text{COOH}$ }





Figure 1: Types of fatty acid



Lipid can be divided according to their chemical composition to

I. Simple lipids: They are esters of fatty acids with various alcohols. They include oils and fats which are esters of fatty acids with glycerol (e.g. **triglycerides**). The triacylglycerol is the simplest and most common fat. It is the form in which lipids are stored in the cell.




Figure 2: A triglyceride is a lipid that is formed by the esterification of glycerol with three fatty acid molecules

II. Complex lipids: Lipids are linking with other compounds, such as phosphate (e.g. **phospholipids**), carbohydrates (e.g. **glycolipids**) and proteins (e.g. **lipoproteins**).

II. Precursor and derived lipids: They include fatty acids, cholesterol, steroid hormones, fat-soluble vitamins (vitamin A, D, E & K) and eicosanoids, which include **prostaglandins, leukotrienes and thromboxanes**.



Biological functions of lipids

1. Lipids provide energy storage and metabolic fuels. They give more energy than carbohydrate and proteins.
2. Cholesterol is a precursor of steroid hormones, bile acids (which help in digestion of dietary fats) and vitamin D.
3. Lipids are found naturally in all living organisms. They have functional and structural components of the cell membrane.
4. Lipoproteins (e.g. LDL & HDL) are a mean for transporting lipids in blood.
5. Imbalance in lipid metabolism can lead to major clinical problems, such as obesity and atherosclerosis.



I. Chemicals reactions of lipids

A- Solubility Test



Objective: to test the solubility of oils in polar & non-polar solvent.

How to perform the test?

- 1- **Add** about **2ml** of the provided oil sample to **5ml** of water in a test tube and **try to mix** oil with water. By shaking, oil and water mix initially; however, they gradually separate out to form 2 layers.
- 2- **Repeat** the experiment using chloroform instead of water.



Result

Oil **dissolves** in the **organic solvent** “chloroform” but **doesn’t dissolve** in **water**, because oils contain long hydrocarbon tails which are hydrophobic in nature.



II- Chemical reaction of fatty acids

Fatty acids: **soluble** in (organic solvents) chloroform and ether, **insoluble** in water

A - Cupper Acetate Test



Objective

To **distinguish** between saturated and unsaturated fatty acid.

How to perform the test?

- 1- **Add 1ml** of each fatty acids (stearic and oleic acids) in 2 test tubes.
- 2- Each of test tube contain **add 3ml** of petroleum ether, shake tubes until fatty acid dissolve.
- 3- **Add 1ml** of cupper acetate ($\text{Cu}_2(\text{CH}_3\text{COO})_4$), **shake** well for minute and wait .



Result



Stearic acid



Oleic acid

The **Blue/green precipitate** is appearing in **the bottom** of tube if fatty acid is **saturated** (palmitic, stearic). But if it is **unsaturated** fatty acids (oleic) the dark **green color** in the **upper solution** will appear.

B – Un saturated Test



Objective

To **indicates** the presence of **double bonds** in the lipid sample.



Principle

Halogens **react** across double bonds .The decolorization of solution of iodine by lipid indicates the presence of double bonds in the lipid sample.

How to perform the test:

- 1- **Add 1ml** of each fatty acids (stearic and oleic acids) in 2 test tubes.
- 2- **Add 1ml** of chloroform, **shake** tubes until fatty acid dissolve.
- 3- **Add drops** iodide alcohol (drop by drop) , **shake** well.

The halogens are non-metallic elements such as fluorine, chlorine, bromine, iodine, astatine



Result

Unsaturated acids (Oleic acid): Brown color of the iodine **disappears**

Saturated acids (Stearic acid): Brown color of the iodine **remains**



III- Chemical reaction of glycerol

Glycerol (or glycerin, glycerine) is a simple **polyol** compound. It is a colourless, odourless, viscous liquid that is widely used as a solvent and in pharmaceutical formulations. Glycerol has three hydrophilic hydroxyl groups (**trihydroxyalcohol**) that are responsible for its **solubility** in water and its hygroscopic nature. The glycerol **backbone** is central to all lipids known as triglycerides. Glycerol is sweet-tasting and of low toxicity.

Polyol (Sugar alcohols):
are organic compounds,
typically derived
from sugars

A- Denestan experiment:

How to perform the test:

- 1- **Put 1ml** of borax solution in clean test tube.
- 2- **Add one drop** of phenolphthalein (ph.ph), the **pink color** is formation.
- 3- **Add** dilute glycerol until the red color disappears.
- 4- **Heat** in the boiling water bath until → **pink color produced again**.

Remember!
b.w.b = boiling
in water bath



Result

Pink color

Glycerol

Disappears
Pink color

b.w.b

Pink color
produced again

Pink color **disappears** and produced **again** after heating

B- Oxidation test



Objective

To **identify** redox reactions and exothermic reactions of glycerol



Principle

The potassium permanganate **oxidizes** the glycerol to **carbon dioxide** and **water** (hence the steam) and is itself reduced.

How to perform the test

- 1- **Put** **1ml** of glycerol solution in test tube.
- 2- **Add** **one ml** of potassium dichromate($K_2Cr_2O_7$).
- 3- **Add** **one drop** of concentrated sulfuric acid (H_2SO_4).
- 4- **Heat** in the boiling water bath few minuets

Potassium permanganate

is a strong oxidizing


agent and **glycerol**

is an easily oxidized

substance



Result



The solution turns to **green color**



VI- Chemical reaction of cholesterol

Cholesterol is a waxy substance, which is made in the body by the liver but is also found in some foods also can be taken in from food. It soluble in chloroform and insoluble in water

A- Salkowski test



Objective

To **detect** the presence of cholesterol



Principle

The cholesterol is react as a typical alcohol with a strong concentrated acids and the product are **colored substances**.

How to perform the test

- 1- **Dissolve** small amount of cholesterol in **2 ml** of chloroform (in a dry test tube)
- 2-**Add 2 ml** of conc. sulfuric acid (H_2SO_4)



Result



The solution becomes **red glow yellow**

B-Liebermann-Burchard test:

Objective

To **detect** the presence of cholesterol


Principle

The cholesterol is react as a typical alcohol with a strong concentrated acids and the product are **colored substances**. Acetic anhydride are used as **solvent** and **dehydrating agents**, and the sulfuric acid is used as **dehydrating** and **oxidizing agent** .A positive result is observed when the solution becomes red , then blue, and finally bluish –green color.


How to perform the test

- 1- **Dissolve** small amount of cholesterol in **2 ml** of chloroform (in a dry test tube)
- 2- **Add 10 drops** of acetic anhydride ($C_4H_6O_3$), mix.
- 3-**Add 2 drops** of conc. sulfuric acid (H_2SO_4)

Result



The solution becomes red , then blue, and finally **bluish–green color**





RESULTS & LAB REPORT

Unknown (fatty acid -stearic acid)

Physical properties

Chemical properties

| Test | Observation | Result |
|------|-------------|--------|
| | | |



RESULTS & LAB REPORT

Unknown (fatty acid -oleic acid)

Physical properties

Chemical properties

| Test | Observation | Result |
|------|-------------|--------|
| | | |



RESULTS & LAB REPORT

Unknown (glycerol)

Physical properties

Chemical properties

| Test | Observation | Result |
|------|-------------|--------|
| | | |



RESULTS & LAB REPORT

Unknown (cholesterol)

Physical properties

Chemical properties

| Test | Observation | Result |
|------|-------------|--------|
| | | |



RESULTS & LAB REPORT

Unknown (1)

Physical properties

Chemical properties

| Test | Observation | Result |
|------|-------------|--------|
| | | |



RESULTS & LAB REPORT

Unknown (2)

Physical properties

Chemical properties

| Test | Observation | Result |
|------|-------------|--------|
| | | |



RESULTS & LAB REPORT

Unknown (3)

Physical properties

Chemical properties

| Test | Observation | Result |
|------|-------------|--------|
| | | |

Qualitative analysis of proteins and amino acids



Introduction of protein

Protein is a complex, high-molecular-weight organic compound that consists of amino acids joined by peptide bonds, its contain C, H, O, and Nitrogen.

Common structure: Central C, with an H, amino group (NH_2), and an acid group (COOH), and a side group (R), proteins made up of about 20 different amino acids. **Unique Side Groups**, differ in size, shape, electrical charge.




Figure 1: The structure of amino acid



The function of protein

Proteins are probably **the most important class** of biochemical molecules, although of course lipids and carbohydrates are also essential for life. They are the **basis** for **the major structural components** of animal and human tissue.

Proteins are **essential** to **the structure** and **function** of all **living cells** and **viruses**. Many proteins are **enzymes** or **subunits of enzymes**, which is catalyzing chemical reactions. Other proteins play structural or mechanical roles, such as those that form the struts and joints of the cytoskeleton, serving as biological scaffolds for the mechanical integrity and tissue signaling functions. It can be **hydrolyzed** by acids, bases or specific enzymes.



The classification of Proteins based on solubility and composition

According to this classification, proteins are divided into three main groups as **simple**, **conjugated** and **derived** proteins.

1- Simple proteins

This group includes proteins **containing only amino acids**, as structural components. These proteins are further classified based on their solubility in different solvents as well as their heat coagulability (The ability to coagulate, of being coagulable). Such as **Egg albumin** and **Globulins**.

2- Conjugated or compound proteins

These are simple proteins **combined** with **some non-protein substances** known as **prosthetic groups**. The nature of the non-protein or prosthetic groups is the basis for the sub classification of conjugated proteins. Such as **nucleoproteins** (combination of simple basic proteins) and **phosphoproteins** (combination with phosphoric acid, e.g. **Casein** of milk).

3- Derived proteins

These are **proteins derived by partial to complete hydrolysis** from the simple or conjugated proteins by the action of **heat, enzymes** or **chemical reagents**. This group also includes the artificially-produced polypeptides. Such as **peptone** and **gelatin**.



Lab experiments.

A- Physical properties

Solubility- color- shape.

Solubility of proteins

Identify the solubility of different types of proteins in different solvents.



Principle

Simple proteins dissolve easily in water while **soluble protein complex** by low solubility in distilled water. The proteins are derived insoluble in water but soluble in alkaline solutions.

How to perform the test:

- 1- Test the solubility of each (**albumin, casein, gelatin, peptone**) in **cold water** and then in **hot water** and then in **sodium hydroxide solution**.
- 2- Record solubility of all proteins in cold water and hot water and sodium hydroxide solution in the following table of results.

| Protein | Albumin | Casein | Gelatin | Peptone |
|-----------------|---------|--------|---------|---------|
| Distilled water | | | | |
| Hot water | | | | |

B. Chemical properties (color test)

1- Biuret Test



Principle

It is the **general test** for all proteins. Biuret reagent is dilute **Copper (II) sulfate** (CuSO_4) in strong alkaline medium. Alkaline CuSO_4 **reacts** with all compounds containing 2 or more peptide bonds to give a **blue-violet color**.

How to perform the test:

1- **Put 2ml** from protein solution (albumin, casein, gelatin, and peptone) in 4 test tubes.

2- **Add 1 ml** of biuret reagent (2 ml of NaOH 10% + 1 drop of CuSO_4 1%).



3- **Mix** well.

Result

Show **positive test** for **all proteins**, is indicated by the **formation of a violet color**.



Positive test



Notes

The **reaction with peptone** may be gave **pink color** when the solution is not concentrated.

Organized by: Duaa Zahim & Dalal Al-Saedi

2- Ninhydrin Test



Principle

Amino acids **contain** a **free amino group** and a **free carboxylic acid group** that **react** together with **ninhydrin** to produce a colored product. When an amino group is attached to the first, or alpha, carbon on the amino acids carbon chain, the amino group's nitrogen atom is part of a **blue-purple product**. Proteins also contain free amino groups on the alphacarbon and can react with ninhydrin to produce a blue-purple product.


How to perform the test:

- 1- **Put 1 ml** from protein solution (albumin, casein, gelatin, and peptone) in 4 test tubes.
- 2- **Add 1 ml** of 0.5% ninhydrin – ethanol reagent solution.
- 3- **Mix well**
- 4- **Place** the test tubes into the boiling-water bath for **5 minute**.



Result

Show positive test for **all proteins**, is indicated by the **formation of a blue-purple**



Positive test



Ninhydrin is a strong oxidizing agent, it should be handled with care, and applied apart from contact with skin or eyes. The gloves and mask is a must be used.

3- Coagulation (clotting) Test



Principle

In this test we see the **denaturing proteins at high temperature** due to change shape, chemical and physical properties of proteins.

How to perform the test

1- **Put 1ml** from protein solution (albumin, casein, gelatin, and peptone) in 4 test tubes.



2- **Heat** the flame directly on the tube and slightly inclined.

Result

Show **positive test** for **albumin only**, is indicated by the formation of **coagulation of the solution**



Positive test

4- Freezing Test



Principle

In this test we see the **change shape at low temperature** due to change shape, chemical and physical properties of proteins.

How to perform the test

- 1- **Put 3ml** from protein solution (albumin, casein, gelatin, and peptone) in 4 test tubes.
- 2- **Put** the tube in the Cup by the snow and leave for **five minutes**.



Result

Show positive test for **gelatin only**, is indicated by the formation of **freezing solution**



Positive test

5- Phosphorus Test



Principle

The basic idea of this test is to **detect** the presence of amino acids that contain the element of **phosphorus**

How to perform the test

- 1- **Put 1ml** from protein solution (albumin, casein, gelatin, and peptone) in 4 test tubes.
- 2 - **Add 1ml** of alkaline solution of 20% NaOH.
- 3 - **Heat** in a boiling- water bath for **10 minutes**.
- 4-**Add** in the order as follows: (**1ml** of Molbidic acid and + **1ml** sodium sulphate + **1ml** of 0.5% hydroquinone) respectively.



Result

Show positive test for casein only, is indicated by the **formation of blue color**.

6- Sulfur Test



Principle

The basic idea of this test is to **detect** the presence of amino acids that contain sulfur such as **cysteine**.

How to perform the test

- 1- **Put 1ml** from protein solution (albumin, casein, gelatin, and peptone) in 4 test tubes.
- 2 - **Add 1ml** of alkaline solution of 40% NaOH.
- 3 - **Heat** in a boiling- water bath for **10 minutes**.
- 4-**Cool** solution and add **1ml** 1% lead acetate.



Result

Show positive test for **albumin**, is indicated by the formation of **black-brown precipitation**.

Show positive test for **peptone**, is indicated by the formation of **yellow precipitation**.

7- Sakaguchi Test



Principle

The basic idea of this test is to **detect** the presence of amino acids that contain **guanidine group** such as **arginine**.

How to perform the test:

- 1- **Put 2ml** from protein solution (albumin, casein, gelatin, and peptone) in 4 test tubes.
- 2- **Add 1ml** of alkaline solution of 10% NaOH.
- 3 - **Add 1ml** of 10% α -naphthol.
- 4 - **Mix** well and cool in ice.
- 5- **Add 1ml** alkaline hypochlorite solution.



Result

Show positive test for **albumin, casein, gelatin and peptone** are indicated by the **formation of color**.

8- Xanthoproteic Test



Principle

The basic idea of this test is to **detect** the presence of amino acids that contain **benzene ring** such **phenylalanine**, **tryptophan** and **tyrosine**.

How to perform the test:

- 1- **Put 2ml** from protein solution (albumin, casein, gelatin, and peptone) in 4 test tubes.
- 2- **Add 1ml** of concentrated nitric acid (HNO_3).
- 3 - **Heat** the mixture in water bath for 30 minutes.
- 4-**Cool** and add drop-wise 40% NaOH to render the solution alkaline.



Result

Show positive test for **albumin, casein and peptone** are indicated by the formation of **orange color or orange precipitation**.



Notes

The gelatin → remains yellow color.

The peptone → convert to orange color.

The albumin and casein → convert to orange precipitation.



RESULTS & LAB REPORT

Unknown (albumin)

Physical properties:

Chemical properties:

| Test | Observation | Result |
|------|-------------|--------|
| | | |



RESULTS & LAB REPORT

Unknown (casein)

Physical properties

Chemical properties

| Test | Observation | Result |
|------|-------------|--------|
| | | |



RESULTS & LAB REPORT

Unknown (gelatin)

Physical properties

Chemical propertie

| Test | Observation | Result |
|------|-------------|--------|
| | | |



RESULTS & LAB REPORT

Unknown (peptone)

Physical properties

Chemical properties

| Test | Observation | Result |
|------|-------------|--------|
| | | |



RESULTS & LAB REPORT

Unknown (1)

Physical properties

Chemical properties

| Test | Observation | Result |
|------|-------------|--------|
| | | |



RESULTS & LAB REPORT

Unknown (2)

Physical properties:

Chemical properties:

| Test | Observation | Result |
|------|-------------|--------|
| | | |



RESULTS & LAB REPORT


Unknown (3)

Physical properties:

Chemical properties:

| Test | Observation | Result |
|------|-------------|--------|
| | | |

Scheme for Unknown Proteins



4 B

Exp: Freezing Test

1 ml of solution, cool in ice



Result: No Freezing

∴ Protein isn't gelatin either

(Casein or Peptone)



Result: Freezing

∴ Protein is Gelatin



5 B

Exp: Phosphorus Test

1 ml of solution add 2 mL of alkaline solution of 20% NaOH. → Δ(b.w.b) for 10 min
Then, add 1 ml of Molbidic acid and + 1 ml sodium sulphate + 1 ml of 0.5% hydroquinone



Result: No Blue color (green)

∴ Protein isn't casein

May be Peptone



Result: Blue color

∴ Protein is Casien

6


Exp: Sulfur Test

1 ml of solution add 1 ml of alkaline solution of 40% NaOH.. → Δ(b.w.b) for 10 min
Then, cool solution and add 1% lead acetate.


Result: Yellow-brown color

∴ Protein is Peptone

Scheme for Unknown Carbohydrates



Scheme for Unknown Lipids





RESULTS & LAB REPORT

Unknown (1)

Physical properties

Chemical properties

| Test | Observation | Result |
|------|-------------|--------|
| | | |



RESULTS & LAB REPORT

Unknown (2)

Physical properties

Chemical properties

| Test | Observation | Result |
|------|-------------|--------|
| | | |



RESULTS & LAB REPORT

Unknown (3)

Physical properties

Chemical properties

| Test | Observation | Result |
|------|-------------|--------|
| | | |